DESCRIPTION

Instant Fortified Wholemilk Powder is a soluble powder made by spray drying fresh pasteurised wholemilk.

FEATURES

- Excellent solubility
- Good flow properties
- Full fat content
- Enriched with Vitamin A and Vitamin D
- Rich creamy flavour
- Contains lecithin

SUGGESTED USES

• Instant Wholemilk Powder is typically re-packed for consumer applications

CERTIFICATION

Halal

ORIGIN

Product of New Zealand



NUTRITIONAL I	NFORMATION
----------------------	------------

TYPICAL

(per 100g of product)

Energy (kJ)	2120
Protein (N x 6.38) (g)	24.0
Moisture (g)	2.8
Fat (g)	28.4
Carbohydrate (g)	39.0
Ash (g)	5.8

Fonterra will only ship this product to countries listed on the Import Eligibility Statement and Bidders Contracting Information file.

INGREDIENTS (ALLERGENS IN BOLD)

Cows' milk, Soy lecithin, Vitamin A, Vitamin D

Country regulations for product labelling vary. Fonterra advises customers that they need to check local regulations to determine the correct labelling of this ingredient.

WHOLEMILK POWDER - Instant Fortified

GDT Specification - Fonterra New Zealand

PRODUCT BULLETIN

PB.560 VERSION:01.0320

UNRESTRICTED

PACKAGING

Multi-wall bag with a paper outer and an inner plastic liner. No staples or metal fasteners are used.

Net Weight: 25 kgGross Weight: 25.35 kg

• Typical Bag Dimension: 95 x 53 x 14 cm



STORAGE AND HANDLING

Wholemilk Powder is hygroscopic and can absorb odours. Therefore adequate protection is essential. It is recommended that product is stored at below 25°C, relative humidity below 65% and in an odour free environment. Stocks should be used in rotation preferably within 24 months of manufacture.



Handle with care



Store in a cool, dry place



TRUSTED QUALITY

COMPLIANCE:

Our manufacturing facilities adhere to the highest global standards in food safety and quality, verified though regular internal and external audit programs and accredited with FSSC22000 certification.

Our products and ingredients are tested throughout the supply chain, from milk collection through to final product release, using internationally recognised test methods to verify quality and compliance.

Our products are securely packed in multi-wall bags with three layers of protection and labelled to enable full trace back, so you can trust in the knowledge it was made with care.



WHOLEMILK POWDER - Instant Fortified

GDT Specification - Fonterra New Zealand

VERSION: 01.0320

NUTRITIONAL ANALYSIS TYPICAL

(p	er 1	00g	of p	prod	uct)

	(per 100g of product)
Energy	2120 kJ
Calories	510 kcal
Fat	28.4 g
Energy from fat	1050 kJ
Calories from fat	250 kcal
Saturated fatty acids	18.8 g
Monounsaturated fatty acids	5.7 g
Polyunsaturated fatty acids	0.4 g
Trans fatty acids	1.3 g
Cholesterol	90 mg
Total Carbohydrate	39.0 g
Lactose	39.0 g
Dietary Fibre	0 g
Calcium	820 mg
Sodium	255 mg
Vitamin A	850 µg
Vitamin D	10 µg



WHOLEMILK POWDER - Instant Fortified GDT Specification - Fonterra New Zealand

PRODUCT SPECIFICATION

PB.560

VERSION: 01.0320

JNRESTRICTED

CHEMICAL	UNITS		MIN		MAX	REFERENCE METHOD
Protein (N x 6.38) as is	% m/m		23.3		26.7	ISO 8968-1
Milk Protein in Milk Solids non fat	% m/m		34			Calculation
Moisture	% m/m				3.5	ISO 5537
Fat	% m/m		28.0		30.0	ISO 1736
Titratable Acidity	% m/v				0.15	SMEDP 15.020
FUNCTIONAL	UNITS		MIN		MAX	REFERENCE METHOD
Wettability [25°C]	seconds				40	ISO/TS 17758
MICROBIOLOGICAL	UNITS		MIN		MAX	REFERENCE METHOD
Aerobic Plate Count	cfu/g				10,000	ISO 4833-1
Coliforms	/g				Not Detected	ISO 11866-1
Escherichia coli	/g				Not Detected	ISO 11866-1
Yeast and Mould	cfu/g				50	ISO 6611
Coagulase Positive Staphylococci	/g				Not Detected	ISO 6888-3
Listeria	/125g				Not Detected	AOAC 051603 /ISO11290
Salmonella	/750g				Not Detected	AOAC 061504 / ISO 6579
PHYSICAL	UNITS		MIN		MAX	REFERENCE METHOD
Bulk Density	g/ml		0.44		0.52	100 Tap
Colour	Typical/Atyp	ical	Typical			Fonterra
Foreign matter	/32.5g		Pass			Fonterra
Scorched particles	/32.5g				В	SMEDP 15.172
Insolubility Index	ml				1.0	ISO 8156
SENSORY	UNITS		MIN		MAX	REFERENCE METHOD
Flavour	Typical/Atypi	cal	Typical			Sensory Analysis
VITAMINS AND MINERALS	UNITS	TYP	ICAL	MIN	MAX	REFERENCE METHOD
Vitamin A	μg/100g			600	1100	HPLC
Vitamin D ₃	μg/100g			7.5	13.0	UPLC-MS/MS
CONTAMINANTS	UNITS		MIN		MAX	REFERENCE METHOD
Inhibitory Substances					Not Detected	DELVO SP-NT

OTHER

Fonterra will only ship this product to countries listed on the Import Eligibility Statement and Bidder Contracting Information File.





WHOLEMILK POWDER, Instant Fortified

GDT Specification - Fonterra New Zealand

IMPORT ELIGIBILITY

DATE: 10 September 2024

UNRESTRICTED

Fonterra will only Ship Instant Fortified Wholemilk Powder GDT Specification – Fonterra New Zealand to the Markets that are specifically listed below:

If you want this product shipped to a country that is not listed below, please contact a Fonterra representative. There is no "obligation" on Fonterra to ship to a country that is not listed or has not been agreed between the customer and Fonterra prior to purchase.

** Prior to contracting with Fonterra importers or agents must ensure that they meet all registration requirements imposed by their importing government agencies.

Destination Country

Bahrain Mongolia Bangladesh Morocco Cambodia Myanmar Cameroon New Caledonia China **Nicaragua Cote d'Ivoire Nigeria Egypt Oman **El Salvador Pakistan **Panama

Georgia Papua New Guinea

Ghana **Philippines** **Guatemala Saudi Arabia Guyana Senegal Singapore Hong Kong Indonesia Solomon Islands South Africa Iraq Jamaica St Vincent Jordan Sudan Kenya Syria Kuwait Taiwan **Thailand Libya

Malaysia United Arab Emirates

Maldives Vanuatu Mauritius Vietnam

Mexico

This Eligibility Statement should be used in conjunction with Product Bulletin Number and Selling Specification Number PB/SS 560

